

Product Datasheet

Human C4orf48 ELISA Kit (orb404926)

Description	Human Neuropeptide-like protein C4orf48(C4orf48) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	C4orf48
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q5BLP8
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.