

Product Datasheet

Human SGSH ELISA Kit (orb404952)

Description	Human N-sulphoglucosamine sulphohydrolase(SGSH) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	SGSH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P51688
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.