

Product Datasheet

Human NOD1 ELISA Kit (orb404964)

Description	Human Nucleotide-binding oligomerization domain-containing protein 1(NOD1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	NOD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y239
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.