

Product Datasheet

Human GPX4 ELISA Kit (orb405019)

Description	Human Phospholipid hydroperoxide glutathione peroxidase, mitochondrial(GPX4) ELISA kit
Reactivity	Human
Range	7.8 µU/mL-500 µU/mL
Target	GPX4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P36969
Sensitivity	1.95 µU/mL
Expiration Date	Please enquire.