

Product Datasheet

Human SERPINB2 ELISA Kit (orb405027)

Description	Human Plasminogen activator inhibitor 2(SERPINB2) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	SERPINB2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P05120
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.