

Product Datasheet

Human C10orf54 ELISA Kit (orb405030)

Description	Human Platelet receptor Gi24(C10orf54) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	C10orf54
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9H7M9
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.