

Product Datasheet

Human PTRF ELISA Kit (orb405040)

Description	Human Polymerase I and transcript release factor(PTRF) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Target	PTRF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6NZI2
Sensitivity	9.4 pg/mL
Expiration Date	Please enquire.