

Product Datasheet

Human WT1 ELISA Kit (orb405052)

Description	Human PRKC apoptosis WT1 regulator protein(PAWR) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	PAWR
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q96IZ0
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.