

Product Datasheet

Human PDCD1LG2 ELISA Kit (orb405069)

Description Human Programmed cell death 1 ligand 2(PDCD1LG2) ELISA kit

Reactivity Human

Range 0.45 ng/ml-30 ng/ml.

Target PDCD1LG2

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates.

Uniprot ID **Q9BQ51**

Sensitivity 0.11 ng/ml

Expiration Date Please enquire.