

## Product Datasheet

### Human EIF2C2 ELISA Kit (orb405094)

|                          |                                                    |
|--------------------------|----------------------------------------------------|
| <b>Description</b>       | Human Protein argonaute-2(EIF2C2) ELISA kit        |
| <b>Reactivity</b>        | Human                                              |
| <b>Range</b>             | 25 pg/mL-1600 pg/mL                                |
| <b>Target</b>            | EIF2C2                                             |
| <b>Note</b>              | For research use only                              |
| <b>Application notes</b> | Detection Wavelength: 450 nm                       |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates, cell lysates |
| <b>Uniprot ID</b>        | <b>Q9UKV8</b>                                      |
| <b>Sensitivity</b>       | 6.25 pg/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                    |