

Product Datasheet

Human NELL1 ELISA Kit (orb405125)

Description	Human Protein kinase C-binding protein NELL1(NELL1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	NELL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q92832
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.