

Product Datasheet

Human WNT5B ELISA Kit (orb405150)

Description	Human Protein Wnt-5b(WNT5B) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	WNT5B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9H1J7
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.