

Product Datasheet

Human TGM1 ELISA Kit (orb405163)

Description	Human Protein-glutamine gamma-glutamyltransferase K(TGM1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	TGM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P22735
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.