

Product Datasheet

Human RIPK3 ELISA Kit (orb405212)

Description	Human Receptor-interacting serine/threonine-protein kinase 3(RIPK3) ELISA kit
Reactivity	Human
Range	0.156 ng/ml-10 ng/ml
Target	RIPK3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Y572
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.