

Product Datasheet

Human PTPRZ1 ELISA Kit (orb405216)

Description	Human Receptor-type tyrosine-protein phosphatase zeta(PTPRZ1) ELISA kit
Reactivity	Human
Range	93.75 pg/mL-6000 pg/mL
Target	PTPRZ1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P23471
Sensitivity	23.43 pg/ml
Expiration Date	Please enquire.