

Product Datasheet

Human RTN4RL2 ELISA Kit (orb405225)

Description	Human Reticulon-4 receptor-like 2(RTN4RL2) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	RTN4RL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q86UN3
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.