

## Product Datasheet

### Human RNF207 ELISA Kit (orb405247)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human RING finger protein 207(RNF207) ELISA kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                           |
| <b>Target</b>            | RNF207  |
| <b>Note</b>              | For research use only                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates |
| <b>Uniprot ID</b>        | <b>Q6ZRF8</b>                                   |
| <b>Sensitivity</b>       | 3.9 pg/mL                                       |
| <b>Expiration Date</b>   | Please enquire.                                 |