

## Product Datasheet

### Human ROBO4 ELISA Kit (orb405252)

<b>Description</b>	Human Roundabout homolog 4(ROBO4) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/mL-4000 pg/mL
<b>Target</b>	ROBO4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernates, urine
<b>Uniprot ID</b>	<a href="#">Q8WZ75</a>
<b>Sensitivity</b>	15.6 pg/mL
<b>Expiration Date</b>	Please enquire.