

Product Datasheet

Human SCGB1D2 ELISA Kit (orb405267)

Description	Human Secretoglobin family 1D member 2(SCGB1D2) ELISA kit
Reactivity	Human
Range	31.25 ng/mL-2000 ng/mL
Target	SCGB1D2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O95969
Sensitivity	7.81 ng/mL
Expiration Date	Please enquire.