

## Product Datasheet

### Human SCG2 ELISA Kit (orb405270)

|                          |                                          |
|--------------------------|------------------------------------------|
| <b>Description</b>       | Human Secretogranin-2(SCG2)<br>ELISA kit |
| <b>Reactivity</b>        | Human                                    |
| <b>Range</b>             | 2.34 ng/ml - 150ng/ml                    |
| <b>Target</b>            | SCG2                                     |
| <b>Note</b>              | For research use only                    |
| <b>Application notes</b> | Detection Wavelength: 450 nm             |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates     |
| <b>Uniprot ID</b>        | <b>P13521</b>                            |
| <b>Sensitivity</b>       | 0.59 ng/ml                               |
| <b>Expiration Date</b>   | Please enquire.                          |