

Product Datasheet

Human SPINK5 ELISA Kit (orb405288)

Description	Human Serine protease inhibitor Kazal-type 5(SPINK5) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	SPINK5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NQ38
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.