

Product Datasheet

Human SFRS1 ELISA Kit (orb405291)

Description	Human Serine/arginine-rich splicing factor 1(SFRS1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	SFRS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q07955
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.