

## Product Datasheet

### Human SIGLEC8 ELISA Kit (orb405322)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human Sialic acid-binding Ig-like lectin 8(SIGLEC8) ELISA kit |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                                          |
| <b>Target</b>            | SIGLEC8                                                       |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates               |
| <b>Uniprot ID</b>        | <b>Q9NYZ4</b>                                                 |
| <b>Sensitivity</b>       | 0.039 ng/mL                                                   |
| <b>Expiration Date</b>   | Please enquire.                                               |