

Product Datasheet

Human SIGLEC9 ELISA Kit (orb405323)

Description Human Sialic acid-binding Ig-like lectin 9(SIGLEC9) ELISA kit

Reactivity Human

Range 125 pg/ml-8000 pg/ml

Target SIGLEC9

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum,plasma,tissue homogenates.

Uniprot ID **Q9Y336**

Sensitivity 31.25 pg/ml

Expiration Date Please enquire.