

Product Datasheet

Human SCUBE1 ELISA Kit (orb405327)

Description	Human Signal Peptide, CUB and EGF-like Domain-containing Protein 1(SCUBE1)ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	SCUBE1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8IWH4
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.