

Product Datasheet

SCAND1 antibody (orb40533)

Description	SCAND1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human SCAN domain-containing protein 1 protein (1-179AA)
Target	SCAND1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57086
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.