

Product Datasheet

Human SPRR2B ELISA Kit (orb405335)

Description	Human Small proline-rich protein 2B(SPRR2B) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	SPRR2B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P35325
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.