

Product Datasheet

Human ATP1A1 ELISA Kit (orb405340)

Description	Human Sodium/potassium-transporting ATPase subunit alpha-1(ATP1A1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	ATP1A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P05023
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.