

Product Datasheet

Human CCT5 ELISA Kit (orb405399)

Description	Human T-complex protein 1 subunit epsilon(CCT5) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	CCT5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates, cell lysates
Uniprot ID	P48643
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.