

## Product Datasheet

### Human CCT6A ELISA Kit (orb405401)

<b>Description</b>	Human T-complex protein 1 subunit zeta(CCT6A) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	25 pg/mL-1600 pg/mL
<b>Target</b>	CCT6A
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P40227</b>
<b>Sensitivity</b>	6.25 pg/mL
<b>Expiration Date</b>	Please enquire.