

Product Datasheet

Human TXNDC5 ELISA Kit (orb405411)

Description	Human Thioredoxin domain-containing protein 5(TXNDC5) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	TXNDC5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8NBS9
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.