

## Product Datasheet

### Human TST ELISA Kit (orb405412)

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Description</b>       | Human Thiosulfate<br>sulfurtransferase(TST) ELISA kit |
| <b>Reactivity</b>        | Human                                                 |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                                 |
| <b>Target</b>            | TST                                                   |
| <b>Note</b>              | For research use only                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell lysates    |
| <b>Uniprot ID</b>        | <b>Q16762</b>                                         |
| <b>Sensitivity</b>       | 3.9 pg/mL                                             |
| <b>Expiration Date</b>   | Please enquire.                                       |