

Product Datasheet

Human THSD4 ELISA Kit (orb405414)

Description	Human Thrombospondin type-1 domain-containing protein 4(THSD4) ELISA kit
Reactivity	Human
Range	12.5 pg/mL-800 pg/mL
Target	THSD4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q6ZMP0
Sensitivity	3.12 pg/mL
Expiration Date	Please enquire.