

Product Datasheet

Human TREM2 ELISA Kit (orb405459)

Description	Human Triggering receptor expressed on myeloid cells 2(TREM2) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	TREM2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NZC2
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.