

Product Datasheet

Human TNFSF14 ELISA Kit (orb405480)

Description	Human Tumor necrosis factor ligand superfamily member 14(TNFSF14) ELISA kit
Reactivity	Human
Range	0.156 ng/ml - 10 ng/ml
Target	TNFSF14
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates.
Uniprot ID	O43557
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.