

Product Datasheet

Human TNFRSF13B ELISA Kit (orb405483)

Description	Human Tumor necrosis factor receptor superfamily member 13B(TNFRSF13B) ELISA kit
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Target	TNFRSF13B
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O14836
Sensitivity	19.5 pg/ml
Expiration Date	Please enquire.