

Product Datasheet

Human TSG101 ELISA Kit (orb405488)

Description	Human Tumor susceptibility gene 101 protein(TSG101) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	TSG101
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q99816
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.