

## Product Datasheet

### Human ABL1 ELISA Kit (orb405491)

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Description</b>       | Human Tyrosine-protein kinase<br>ABL1(ABL1) ELISA kit |
| <b>Reactivity</b>        | Human                                                 |
| <b>Range</b>             | 25 pg/mL-1600 pg/mL                                   |
| <b>Target</b>            | ABL1                                                  |
| <b>Note</b>              | For research use only                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                          |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                     |
| <b>Uniprot ID</b>        | <b>P00519</b>                                         |
| <b>Sensitivity</b>       | 6.25 pg/mL                                            |
| <b>Expiration Date</b>   | Please enquire.                                       |