




Product Datasheet

FCRL4 antibody (orb40551)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	FCRL4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57224
Uniprot ID	Q96PJ5
NCBI	NP_112572.1
Dilution Range	WB: 1:2000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of FCRL4 antibody...



Immunohistochemical analysis of formalin...