

Product Datasheet

Human ATP6V0A2 ELISA Kit (orb405548)

Description	Human V-type proton ATPase 116 kDa subunit a isoform 2(ATP6V0A2) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	ATP6V0A2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Y487
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.