

Product Datasheet

Human XIRP1 ELISA Kit (orb405559)

Description	Human Xin actin-binding repeat-containing protein 1(XIRP1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	XIRP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q702N8
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.