

Product Datasheet

Human CYP24A1 ELISA Kit (orb405577)

Description	Human 1,25-dihydroxyvitamin D(3) 24-hydroxylase, mitochondrial (CYP24A1) ELISA kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	CYP24A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q07973
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.