

Product Datasheet

Human YWHAB ELISA Kit (orb405579)

Description	Human 14-3-3 protein beta/alpha (YWHAB) ELISA kit
Reactivity	Human
Range	6.25 ng/mL-400 ng/mL
Target	YWHAB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P31946
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.