

## Product Datasheet

### Human YWHAH ELISA Kit (orb405580)

|                          |                                           |
|--------------------------|-------------------------------------------|
| <b>Description</b>       | Human 14-3-3 protein eta(YWHAH) ELISA kit |
| <b>Reactivity</b>        | Human                                     |
| <b>Range</b>             | 0.625 ng/mL-40 ng/mL                      |
| <b>Target</b>            | YWHAH                                     |
| <b>Note</b>              | For research use only                     |
| <b>Application notes</b> | Detection Wavelength: 450 nm              |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates         |
| <b>Uniprot ID</b>        | <b>Q04917</b>                             |
| <b>Sensitivity</b>       | 0.156 ng/mL                               |
| <b>Expiration Date</b>   | Please enquire.                           |