

Product Datasheet

Human CYP27B1 ELISA Kit (orb405588)

Description	Human 25-hydroxyvitamin D-1 alpha hydroxylase, mitochondrial(CYP27B1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	CYP27B1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O15528
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.