

Product Datasheet

Human SRD5A1 ELISA Kit (orb405590)

Description	Human 3-oxo-5-alpha-steroid 4-dehydrogenase 1(SRD5A1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	SRD5A1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P18405
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.