

Product Datasheet

Human SRD5A2 ELISA Kit (orb405593)

Description	Human 3-oxo-5-alpha-steroid 4-dehydrogenase 2(SRD5A2) ELISA kit
Reactivity	Human
Range	12 pg/mL-2400 pg/mL
Target	SRD5A2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P31213
Sensitivity	6 pg/mL
Expiration Date	Please enquire.