

Product Datasheet

Human 8-epi-PGF2 α ELISA Kit (orb405604)

Description	Human 8-epi-prostaglandin F2alpha,8-epi-PGF2 α ELISA Kit
Reactivity	Human
Range	78 pg/mL-5000 pg/mL
Target	8-epi-PGF2 α
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, urine, tissue homogenates
Sensitivity	19.5 pg/mL
Expiration Date	Please enquire.