

## Product Datasheet

### Human ADAMTS4 ELISA Kit (orb405624)

|                          |                                                                                 |
|--------------------------|---------------------------------------------------------------------------------|
| <b>Description</b>       | human ADAM metallopeptidase with thrombospondin type 1 motif4,ADAMTS4 ELISA Kit |
| <b>Reactivity</b>        | Human                                                                           |
| <b>Range</b>             | 12.5 ng/ml - 800 ng/ml                                                          |
| <b>Target</b>            | ADAMTS4                                                                         |
| <b>Note</b>              | For research use only                                                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                               |
| <b>Uniprot ID</b>        | <b>075173</b>                                                                   |
| <b>Sensitivity</b>       | 3.125 ng/ml                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                 |