

Product Datasheet

Human AKR1B10 ELISA Kit (orb405626)

Description	Human Aldo-keto reductase family 1 member B10(AKR1B10) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	AKR1B10
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O60218
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.